How to perform event building in a MINIBALL experiment

R. Lutter

May 20, 2005

Abstract

This document describes how to perform event building and to generate ${\tt ROOT}$ trees from med data in a Miniball experiment.

Contents

1	Online data acquisition	2
2	Filling hitBuffers from raw data 2.1 Extracting hits from DGF-4C data	
3	Event building 3.1 Inserting time stamps in ADC/TDC data	4 5 6
4	Replay session	8

1 Online data acquisition

Any data in a MINIBALL experiment will be produced by connecting to the MBS subsystem and reading raw data from the harware modules in VME and CAMAC crates [1]. Online data are then written to disk using the MED format (MBS event data) [2].

Fig. 1 shows data flow as well as software components used to acquire experimental data.

Following files are needed for this step:

• DgfReadout.c user's readout function given by hardware definitions in Config.C auto-generated during config step

- DgfAnalyze.cxx methods and functions to setup the ROOT system and to decode raw data from MBS, auto-generated during the config step
- udef/ProcessEvent.cxx consistency checks and accumulation of diagnostic histograms, has to be provided by user




Figure 1: Online data acquisition

- 2 Filling hitBuffers from raw data
- 2.1 Extracting hits from DGF-4C data
- 2.2 Extracting hits from CAEN ADC/TDC data

3 Event building

In a second step one may create meaningful events from raw data read from MED files. These events may then be written to ROOT trees and used in replay sessions afterwards. Fig. 2 shows how to do the event building.




Figure 2: Event building

Steps to be performed for the event building process:

- in SetCppIfdefs.C : activate event building

 dgf->MakeDefined("__EVENT_BUILDING_ON__", kTRUE, "Start event building");
- in Config.C: define an event class to be used for event building

dgf->IncludeUserClass("udef", "MbEvent.cxx", kTRUE);

when doing the config step for the first time this will generate a template class MbEvent; you then have to modify this class according to your needs

- provide an event building algorithm: TUsrEvtReadout::BuildEvent() in udef/BuildEvent.cxx (table 2, fig. 3)
- in .rootrc file: add an entry TMrbAnalyze.CoincWindow: N
- perform config, compile and link the whole stuff, run C_analyze

3.1 Inserting time stamps in ADC/TDC data

```
Bool_t TUsrEvtReadout::InsertTimeStamp(Int_t TsIdx, Int_t TsChannel, Int_t StartIdx, Int_t NofHbx) {
                                                                  [C++ METHOD]
// Name:
// Purpose:
                 TUsrEvtReadout::InsertTimeStamp
                 Insert time stamps into adc data
                 Int_t TsIdx -- index of timestamping buffer
Int_t TsChannel -- number of timestamping channel
Int_t StartIdx -- index of adc buffer to start with
Int_t NofHbx -- number of buffers to process
// Arguments:
                 kTRUE/kFALSE
// Description:
                 Inserts time stamps from timestamping dgf into adc buffers.
                 Time stamps in dgf buffer are sorted by time,
                 adc/tdc buffers contain single channel data grouped by event number.
                 This algorithm inserts time stamps one by one in adc channel data;
                 channels belonging to one event (i.e. having the same event number)
                 will be marked with same time stamp.
// Keywords:
TUsrHBX *tsHbx = this->GetHBX(TsIdx);
                                                          // connect to hitbuffer of time stamping dgf
 for (Int_t i = StartIdx; i < StartIdx + NofHbx; i++) { // loop over adcs/tdcs</pre>
    tsHbx->ResetIndex();
                                                          // reset to first item in ts-dgf buffer
                                                          // fetch first item (= time stamp)
    TUsrHit * tsHit = tsHbx->FindHit(TsChannel);
    if(!tsHit) break;
                                                          // no time stamps in this buffer (should never be)
    TUsrHBX * cHbx = this->GetHBX(i);
                                                          // connect to current adc/tdc buffer
    cHbx->ResetIndex();
                                                          // reset to first item in buffer
    TUsrHit * cHit = cHbx->NextHit();
                                                          // fetch first adc/tdc hit (= channel data)
    while(tsHit && cHit) {
                                                          // step thru buffers
                                                          // as long as event number doesn't change:
// insert current time
      Int_t evtNo = cHit->GetEventNumber();
                                                          // save current event number from adc/tdc hit
      while(cHit->GetEventNumber() == evtNo) {
       cHit->SetChannelTime(tsHit);
       cHit = cHbx->NextHit();
       if (!cHit) break;
                                                          // end of adc/tdc buffer reached
      tsHit = tsHbx->FindHit(TsChannel);
                                                          // event number in adc/tdc buffer has changed:
                                                                get next time stamp from ts-dgf buffer
 return(kTRUE);
```

Table 1: Algorithm used for time stamp insertion

3.2 Event building algorithm

```
Bool_t TUsrEvtReadout::BuildEvent() {
                                                              [C++ METHOD]
..
.......
               TUsrEvtReadout::BuildEvent
// Name:
// Purpose:
               Event building
// Arguments:
// Results:
               kTRUE/kFALSE
// Exceptions:
// Description: Loops over all hit buffers,
              collects hits within window 'coincWindow'
               stores hits in events of type MbEvent
              then calls method MbEvent::WriteToTree() for each event
              Strategy:
               a binary tree TBtree is used to compare events
                method TBtree::FindObject() calls MbEvent::Compare()
                MbEvent::IsEqual() then tests if time stamps are 'equal' within 'coincWindow'
                hits belonging to an event are stored in a hit buffer of type TClonesArray,
                to get the benefits of TClonesArray these hit buffers have to be static
                they will be allocated only once and reused afterwards.
TBtree * evts = &evtsS;
 TObjArrayIter epI(&poolOfMbEvents);
  while(MbEvent * evt = (MbEvent *) epI.Next()) evt->Reset();
// fill binary tree, create events if necessary
  Int_t evtCount = 0;
 for (Int_t i = kIdxCluster1; i < kIdxCluster1 + kNofSevts; i++) { // loop over all hit buffers
   if (hbx) {
     Int_t hitNo = 0;
     tmpEvent.SetEventTime(hit->GetChannelTime());
// time stamp = 0 -> junk data
// insert time
                                                          // insert time stamp into temp event
       MbEvent * evt = (MbEvent *) evts->FindObject(&tmpEvent); // compare it with events already there if (evt == NULL) { // no match if (evtCount < poolOfMbEvents.GetEntriesFast()) { // try to get a event from pool
           evt = (MbEvent *) poolOfMbEvents[evtCount];
         } else {
           evt = new MbEvent();
                                             // create a new one
                                            // add it to pool
           poolOfMbEvents.Add(evt);
         evt->Reset():
                                             // reset event time, hit buffer etc
         evt->SetEventTime(hit->GetChannelTime()); // insert time stamp
evt->SetCoincWindow(coincWdw); // should know about time window
                                            // add new event to binary tree
         evts->Add(evt);
         evtCount++;
       evt->AddHit(hit):
                                           // append current hit to event
       hitNo++;
// call method MbEvent::WriteToTree() for all events
TBtreeIter evI(evts);
while(MbEvent * evt = (MbEvent *) evI.Next()) evt->WriteToTree(); // write event data to tree
evtsS.Clear();
                                            // clear entries in binary tree
return(kTRUE);
```

Table 2: Algorithm used for the event building process




Figure 3: Event builder: schematic diagram

4 Replay session

Once events have been built and tree data have been written to ROOT file one may start a replay session by reading these pre-sorted data. Fig. 4 shows the data flow for such a replay session. The gain in speed is about a factor of 10 as compared to replaying raw data from MED files.




Figure 4: Offline replay session

User has to provide a method MbEvent::Analyze() containing his analysis. A wrapper class MbExp may be used to establish pointer arrays pointing to different event components such as cores, segments, channels, etc (table 3).

```
[C++ CLASS DEFINITION]
// Name:
                MbExp
// Purpose:
                Wrapper class to access MbEvent data
// Description:
                Organizes MbEvent data in a more "experiment-oriented" way.
                Data may be accessed via cluster, core, and segment numbers (dgf),
                via adc and channel numbers, resp.
                Performs boundary check if requested.
                Indices start with 0.
// Keywords:
class MbExp {
 public:
   MbExp(MbEvent * Event = NULL) { Init(Event); };
                                               // default ctor
   virtual ~MbExp() {};
                                   // default dtor
   TUsrHit * Core(Int_t Cluster, Int_t Core, Bool_t BoundaryCheck = kFALSE);
   TUsrHit * Segment(Int_t Cluster, Int_t Core, Int_t Seg, Bool_t BoundaryCheck = kFALSE);
   TUsrHit * Particle(Int_t Dgf, Int_t Channel, Bool_t BoundaryCheck = kFALSE);
   TUsrHit * Time(Int_t Tdc, Int_t Channel, Bool_t BoundaryCheck = kFALSE);
   void Init(MbEvent * Event = NULL); // initialize ptr arrays with data from MbEvent
 protected:
   TUsrHit * fCores[kNofDgfClusters][kNofDgfCoresPerCluster];
   TUsrHit * fSegs[kNofDgfClusters][kNofDgfCoresPerCluster][kNofDgfSegsPerCore];
   TUsrHit * fPart[kNofAdcs][kNofChannelsPerAdc];
   TUsrHit * fTime[kNofTdcs][kNofChannelsPerTdc];
 protected:
   TUsrHit fEmptyHit;
```

Table 3: Wrapper class to access event components

References

- [1] See also: "Instructions how to use the DAQ in a MINIBALL experiment", http://www.bl.physik.uni-muenchen.de/marabou/html/doc/DaqInstructions.pdf
- [2] For details see: "MED Data Structure", http://www.bl.physik.uni-muenchen.de/marabou/html/doc/MedStructure.pdf